








Required cables for connecting various defibrillators to our live defibrillator compatible manikins:

- Resusci Anne CPR-D (32002xxx, 32005xxx)
- Resusci Anne CPR-D w/airway head (199-00050, 199-05050)
- Resusci Anne Simulator (150-000xx)
- Resusci Anne SkillTrainer (151-000xx)
- ALS SkillMaster (290012xx)
- SimMan (211-00050)
- SimBaby (245-05050)
- ALS Simulator (205-05050)
- MegaCode Kelly (200-05050, 203-05050)
- MegaCode Kid (231-05050, 231-15050)

Defibrillator type	Adapter	Training cable
Heartstart MRx Heartstart XL, 4000* Heartstart FR, FR2, FR2+, FRx	No	Cat no. 945004
Zoll 1200, 1400, 1600 and M-series defibrillators -or other types using similar electrodes	Laerdal Cat no. 05-10100 	
Medtronic LP12 -or other types using Quick Combo connector	Laerdal Cat no. 05-10000 	

<i>*HeartStart MRx, 4000, XL and Codemaster 100, which is equipped with the white barrel-style connector, must use the Heartstart Pads Adapter 05-10200 in addition to the training cable.</i>	Laerdal Cat no. 05-10200 	Cat no. 945004 
<i>HeartStart 3000 needs no training cables.</i>	No	No